

Daniels Accessories

Maximising the efficiency
of your Daniels solutions



MAKING HEALTHCARE SAFER
SOINS DE SANTÉ PLUS SÉCURITAIRES

Position. Storage. **Movement**

Engineered to optimise space, overcome logistical challenges, and enable point-of-use disposal, the Daniels suite of clinical products and accessories are designed as an end-to-end solution, providing safe and seamless transitions between different work environments and applications.

The three signature designs exclusive to Daniels Health which enable logistical versatility and space optimisation are the Daniels peg system, mounting bracket and rail locks. These three systems form the genesis of all our accessory designs.

Daniels Peg System




Mounting Bracket





Rail Locks



Our interchangeable SharpSmart & Accessory systems are designed not only to achieve optimum versatility but also to champion safe handling for hospital staff and reduce the risk of infection. All of our mounting and logistical accessories ensure that SharpSmart collectors are always suspended from the ground, and the maneuverability of our transport accessories take the risk out of waste handling.











The Daniels SharpSmart collector are not just a sharps waste receptacle, they are part of an engineered risk-averse solution that manages waste from cradle to grave with the highest safety and environmental standards ever evidenced in the medical industry.



Sharpsmart Mounting Options

The Daniels suite of bracketing and transportation accessories have been engineered to achieve maximum versatility, optimise space, and overcome logistical challenges within clinical environments. Designed with a small footprint to fit into space-restrictive areas, the Daniels trolleys and floorstands are fitted with a standard mounting bracket which gives optimum flexibility for point-of-use disposal and efficient movement of collectors.

ORDER CODE

7111	STANDARD CRADLE Easy to install design simply hooks over rail on trolley to allow Sharpsmart collectors to face inwards or outwards. Rugged, non-rust stainless steel construction. Compatibility: S14, S22 Regular and Access Plus sizes	
7110	MOUNTED CRADLE Identical concept and design as Standard Cradle, but with additional mounting plate to attach cradle to a wide variety of trolley designs. Compatibility: S14, S22 Regular and Access Plus sizes	
7115	STANDARD WALL BRACKET Universal bracket enabling wall mounting of Sharpsmart range. Bracket kit comes with screws and hole-plugs. For wall mounting of Access Plus Sharpsmarts, see Access Plus bracket. Compatibility: S14, S22 and S32 Regular sizes	
7105	WALL BRACKET - ACCESS PLUS Designed for wall mounting of the Sharpsmart Access Plus range to allow 90-degree lid opening. Compatibility: S14, S22, S32 Access Plus sizes	
7112 Left Side Lock	SHARPSMART LOCKING WALL BRACKET For use in areas where secure, key lock mounting of Sharpsmart collectors is required.	
7113 Right Side Lock	Compatibility: S14, S22 and S32 Regular sizes	
7114	LOCKING MOUNTSMART Designed for versatile inward or outward mounting of Sharpsmart collectors to allow 90-degree lid opening. Compatibility: S14, S22, S32 Regular and Access Plus sizes	
7116	BENCH BRACKET FOR SHARPSMART S14 Provides secure stabilization of S14 collectors on a benchtop. Compatibility: S14 Sharpsmart Regular and Access Plus sizes	
7116	BENCH BRACKET FOR SHARPSMART S22 Provides secure stabilization of S22 collectors on a benchtop. Compatibility: S22 Sharpsmart Regular and Access Plus sizes	
7320	S-SERIES FLOOR STAND Designed for stable floor mounting. Stainless steel and plastic construction is tough and easy to clean. Compatibility: S14, S22 and S32 Regular and Access Plus sizes	
7330	CARTSMART Enables convenient, mobile, off the floor mounting of Sharpsmarts. Plastic construction is lightweight and easy to clean. Compatibility: All Sharpsmart Regular and Access Plus sizes	

Sharpsmart Storage & Logistics

6005019001 MRI CARTSMART

Designed for use in MRI suites and any applications which require non-magnetised equipment. This product has an MR Conditional Rating.

Compatibility: All Sharpsmart Regular and Access Plus sizes



7352 SOFT CLOSE ACCESSMART

Designed for optimized ward movement with an odour-reduction mechanism, this foot-pedal operated trolley can also be fitted with a side mounting bracket to hold Sharpsmart S14 and S22 regular collectors.

Accessmart Trolley Compatibility: 64 Series collectors

S-Series Side Bracket Compatibility: S14, S22 regular only



4004806000 WALLSMART STORAGE SYSTEM

750mm wide

Unique hook system that enables Sharpsmart collectors to be wall mounted for efficient storage within hospitals. Modular system of panels can be arranged to suit different sized walls and storage areas.

7800

1120mm wide

Panel size: 70.9" (H) | Two width sizes are available: 29.5" & 44.1"

Compatibility: All Sharpsmart Regular and Access Plus sizes

The 750mm wide can hold: 15 x S14; or 12 x S22; or 9 x S32; or 4 x 64-Series collectors, or a combination of various sizes.
The 1120mm wide can hold: 25 x S14, or 20 x S22; or 15 x S32; or 6 x 64-Series collectors, or a combination of various sizes.



4017700000 WALLSMART MINI

A wall mounted system for holding up to 3 Sharpsmart collectors.

Standard panel size: 6.0" (H) x 26.5" (W) x 1.5"(D)

Compatibility: All Sharpsmart Regular and Access Plus sizes



6001606002 SMALL INTERNAL DELIVERY CART

Custom made delivery cart with rubber buffers and brakes for movement of Sharpsmarts within a healthcare facility. Holds 40 x S14; or 20 x S22 plus 10 x S14; or 20 x S32, or a combination of various sizes.

S-Series Cart Compatibility: S14, S22, S32 Regular and Access Plus Sizes



7379 SMALL 64 INTERNAL DELIVERY CART

Added compatibility with 64 series collectors only.

7382 LARGE INTERNAL DELIVERY CART

Custom made delivery cart with rubber buffers and brakes for movement of Sharpsmarts within a healthcare facility. Holds 80 x S14; or 48 x S22 plus 8 x S14; or 32 x S32 plus 8 x S22; or 16 x S14, or a combination of various sizes.

S-Series Cart Compatibility: S14, S22, S32 Regular and Access Plus Sizes



7380 LARGE 64 SERIES INTERNAL DELIVERY CART

Added compatibility with 64 series collectors.

4001705000 ADAPTER BRACKET

The 64 series adapter bracket is designed to mount onto the horizontal bars of the delivery cart to facilitate transport of the S64 (and other 64 series) collectors. Each bracket mounts two S64

Compatibility: 64 Series Accessplus only





MAKING HEALTHCARE SAFER
SOINS DE SANTÉ PLUS SÉCURITAIRES

Daniels has operations in the Canada, USA, Australia, New Zealand, UK, Europe and South Africa and is widely regarded as a leader of innovation in medical waste management. Through education, innovation and service, we are committed to making the healthcare industry safer for the people who work in it, and continually driving the development of ecologically sustainable solutions to reduce the environmental impact of healthcare generated waste.